Cleaning Tinkerine Ear Relief

Step 1. Precautions for cleaning

- Before you begin, protect yourself by wearing gloves, a face mask, and a face shield if possible.
- Clean and disinfect the working surface thoroughly.
- Make sure Tinkerine Ear Relief is detached from the face mask.
- Only clean one Tinkerine Ear Relief at a time to avoid cross-contamination.

Step 2. Cleaning the ear relief

- It is recommended to fully submerge Tinkerine Ear Relief in a disinfectant solution to ensure complete sanitization before and after every use.
- After removing Tinkerine Ear Relief from the solution, allow a few minutes for the ear relief to dry. Wipe down any remaining disinfectant droplets with a clean paper towel.
- If submerging is impossible, thoroughly spraying the ear relief with a disinfectant can also be effective. Allow a few minutes to dry. Wipe down any remaining droplets.
- Our recommended disinfectants are:
 - Isopropyl alcohol (70%)
 - Hydrogen peroxide (28%)
- To find out more about specific cleaning products, please refer to the [list of disinfectants approved by Health Canada](#).

Tinkerine Ear Relief is not sterile. Since the conditions and methods of assembly and use of Tinkerine Ear Relief are beyond the manufacturer's control, the manufacturer expressly disclaim any and all liability arising from any assembly and use of Tinkerine Ear Relief. NO WARRANTY OF FITNESS FOR ANY PARTICULAR PURPOSE, WARRANTY OF MERCHANTABILITY OR ANY OTHER WARRANTY, EXPRESSED OR IMPLIED, IS MADE CONCERNING TINKERINE EAR RELIEF. The user should thoroughly test any application before use. By assembling and using Tinkerine Ear Relief, the user does so at their own risk and waives all claims against the manufacturer. The user releases the manufacturer from all liability, whether due to foreseeable or unforeseeable risks, relating to Tinkerine Ear Relief. The user is solely responsible for their assembly and use of the Tinkerine Ear Relief. The manufacturer is not responsible for errors or omissions relating to Tinkerine Ear Relief. The manufacturer does not provide any medical advice.

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